



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High Performance Static Dissipative Cleanroom Sleeves

**Product Code :**  
LBNY-00017-11000265



**Description :**

Features elastic trimmed ends and tapered length for proper fit  
Suitable for exposure to bloodborne pathogens, non-hazardous liquids and any wet general maintenance areas  
Non-linting Environstar-MP fabric is constructed from unique laminated and breathable membrane  
Lint free polyethylene and polypropylene sleeves for use in critical, ESD and low particulate clean environments.  
Length: 18 in.  
White in color.

**Technical Specification :**

High Performance Static Dissipative Cleanroom Sleeves



**Laboratory Instrument India**