



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Safety Apparel Cryogenic Apron

**Product Code :**  
LBNY-00017-11000266



**Description :**

Constructed of Teflon Laminated Nylon  
National Safety Apparel Cryogenic Apron provides waterproof, thermal protection against liquid nitrogen splashes and other extremely cold temperature tasks.  
Water resistant  
Adjustable strap and buckle closure at waist and neck  
Scotchlite Thinsulate liner.

**Technical Specification :**

Safety Apparel Cryogenic Apron



**Laboratory Instrument India**