



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Apparel Cryogenic Apron

Product Code :
LBNY-00017-11000266



Description :

Constructed of Teflon Laminated Nylon
National Safety Apparel Cryogenic Apron provides waterproof, thermal protection against liquid nitrogen splashes and other extremely cold temperature tasks.
Water resistant
Adjustable strap and buckle closure at waist and neck
Scotchlite Thinsulate liner.

Technical Specification :

Safety Apparel Cryogenic Apron



Laboratory Instrument India