



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
IRIS Shoe Covers

**Product Code :**  
LBNY-00017-11000268



**Description :**

Low linting, conforms to cleanrooms up to level ISO 5  
Seamless sole and sides  
Disposable  
Strong and durable  
Waterproof  
Anti-skid  
Disposable shoe covers are made of a proprietary fabric that offers a balance between strength, anti-skid, fluid resistance, and wear resistance.  
Barrier protection  
Color: Light Blue  
Available in universal and XL Sizes.

**Technical Specification :**

IRIS Shoe Covers



**Laboratory Instrument India**

