



Email : sales@laboratoryinstrumentindia.com

Product Name :
IRIS Shoe Covers

Product Code :
LBNY-00017-11000268



Description :

Low linting, conforms to cleanrooms up to level ISO 5
Seamless sole and sides
Disposable
Strong and durable
Waterproof
Anti-skid
Disposable shoe covers are made of a proprietary fabric that offers a balance between strength, anti-skid, fluid resistance, and wear resistance.
Barrier protection
Color: Light Blue
Available in universal and XL Sizes.

Technical Specification :

IRIS Shoe Covers



Laboratory Instrument India

