



Email : sales@laboratoryinstrumentindia.com

Product Name :
Shoe Covers with Non-Skid Print

Product Code :
LBNY-00017-11000270



Description :

Ideal use in medical and pharmaceutical environments, GMP facilities, pre-cleanroom gowning zones and general hygiene

High quality fabric and strong welding construction

Disposable non-skid shoe covers are lightweight, breathable and provide splash protection. Ideal for transitioning between clean and unclean areas. Made of a polypropylene fabric

Non-Sterile, Non-Skid

Shoe covers are 100% latex and silicone free

Choose from a variety of sizes and packaging options.

Triple-elastic runs that allow for a smooth wrap around the foot to prevent tripping

Available in 2 sizes: Universal to fit most shoe sizes, XL for bulkier or steel toe boots

Available in several packaging options bags or dispensing boxes.

Technical Specification :

Shoe Covers with Non-Skid Print



Laboratory Instrument India

