



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Reusable Barrier Lab Coat

Product Code :
LBNY-00017-11000272



Description :

Quick-release snap front closure
Two front pockets and one left breast pocket
Side seam openings for pants pocket access
These specialized lab coats offer excellent fluid resistance while maintaining their breathability.
Splash Resistant Polyester front and sleeves
Breathable Poly Cotton back panel for ventilation
100% polyester fluid-resistant front and sleeves
Industrial launderable
Can be used as part of an Compliance Plan.

Technical Specification :

Protective Wear Reusable Barrier Lab Coat



Laboratory Instrument India