



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Flame-Resistant Lab Coat

Product Code :
LBNY-00017-11000273



Description :

Quick-release snap front
Safety-sleeve cuff with tunneled elastic
Breast pocket with pencil slot
360° fire-resistant identification patches
2 oversized hip pockets
Vent openings on the sides allow access to pants pockets
Flame-Resistant Lab Coats are durable enough to withstand industrial laundering while maintaining fire-resistant properties that are guaranteed for the life of the garment.
Fire-resistant properties guaranteed for life of garment.

Technical Specification :

Protective Wear Flame-Resistant Lab Coat



Laboratory Instrument India