



Email : sales@laboratoryinstrumentindia.com

Product Name :
Aseptic Extra-Safe Lab Jackets — Cranberry

Product Code :
LBNY-00017-11000275



Description :

Autoclavable and washable up to five times
Latex- free, fire retardant
Static-free — no clinging
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.
Hip length
Fluid- resistant, breathable
Snap front, three pockets, knit cuffs and collar.

Technical Specification :

Aseptic Extra-Safe Lab Jackets — Cranberry



Laboratory Instrument India