



Email : sales@laboratoryinstrumentindia.com

Product Name :
Functional Distillation Column

Product Code :
LBNY-0001-990009



Description :

Functional Distillation Column

Technical Specification :

Fractional distillation is particularly suitable for critical applications that require separation of complex mixtures such as fine chemicals, pharmaceuticals, foods, flavors, cosmetics, extracts, etc. – anywhere where a finer level of purity is needed or when component's boiling points are close.

Fractional distillation equipment is available in a multiple range of sizes from the lab scale solvent recovery through processing scale. Pope specializes in mid-size fractionation columns.

These precision designed fractionation units fill the gap between small laboratory bench scale fractional distillation units and large process systems for critical applications. The 1 – 100 gal/hr throughput capacity is ideal for pilot plants and for applications requiring small runs.



Laboratory Instrument India