



Email : sales@laboratoryinstrumentindia.com

Product Name :
Lab Coats — Sky Blue

Product Code :
LBNY-00017-11000276



Description :

Fluid-resistant, breathable
Snap front, three pockets, knit cuffs and collar
Latex-free, fire retardant
Static-free — no clinging
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.
Knee length
Autoclavable and washable up to five times.

Technical Specification :

Lab Coats — Sky Blue



Laboratory Instrument India