



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Protective Wear Reusable Barrier Lab Coat

**Product Code :**  
LBNY-00017-11000278



**Description :**

These specialized lab coats offer excellent fluid resistance while maintaining their breathability.  
100% polyester fluid-resistant front and sleeves  
Quick-release snap front closure  
Breathable Poly Cotton back panel for ventilation  
Splash Resistant Polyester front and sleeves  
Industrial launderable  
Two front pockets and one left breast pocket  
Side seam openings for pants pocket access.

**Technical Specification :**

Protective Wear Reusable Barrier Lab Coat



**Laboratory Instrument India**