



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Reusable Barrier Lab Coat

Product Code :
LBNY-00017-11000278



Description :

These specialized lab coats offer excellent fluid resistance while maintaining their breathability.
100% polyester fluid-resistant front and sleeves
Quick-release snap front closure
Breathable Poly Cotton back panel for ventilation
Splash Resistant Polyester front and sleeves
Industrial launderable
Two front pockets and one left breast pocket
Side seam openings for pants pocket access.

Technical Specification :

Protective Wear Reusable Barrier Lab Coat



Laboratory Instrument India