



Email : sales@laboratoryinstrumentindia.com

Product Name :
Embossed Disposable Gowns

Product Code :
LBNY-00017-11000279



Description :

Large blue gowns feature easy tear-away for quick gown removal and safe disposal
Easy tear-away feature for quick gown removal and safe disposal
Contoured thumb loops secure gown sleeve for added protection
Redirects stress to the film material not the thermal seals to minimize breakage potential
Full-length open back affords better comfort and shoulder coverage
Provides liquid and aerosol blood resistance
Lightweight, non-linting construction
In compliance with FDA requirements and tested in accordance with ASTM standards including breaking load, tensile strength and elongation
Successfully tested for liquid and aerosol blood resistance under ASTM F 1670 In compliance with bloodborne pathogen standards levels 1 & 2 barrier protection
Superior physical properties over imports
Synthetic Blood Penetration Resistance and ASTM F 1671 Bacteriophage Penetration Resistance.

Technical Specification :

Embossed Disposable Gowns

Laboratory Instrument India

