



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Embossed Disposable Gowns

**Product Code :**  
LBNY-00017-11000279



**Description :**

Large blue gowns feature easy tear-away for quick gown removal and safe disposal  
Easy tear-away feature for quick gown removal and safe disposal  
Contoured thumb loops secure gown sleeve for added protection  
Redirects stress to the film material not the thermal seals to minimize breakage potential  
Full-length open back affords better comfort and shoulder coverage  
Provides liquid and aerosol blood resistance  
Lightweight, non-linting construction  
In compliance with FDA requirements and tested in accordance with ASTM standards including breaking load, tensile strength and elongation  
Successfully tested for liquid and aerosol blood resistance under ASTM F 1670 In compliance with bloodborne pathogen standards levels 1 & 2 barrier protection  
Superior physical properties over imports  
Synthetic Blood Penetration Resistance and ASTM F 1671 Bacteriophage Penetration Resistance.

**Technical Specification :**

Embossed Disposable Gowns

---

# Laboratory Instrument India

