



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Blue Premium Bouffant Cap

**Product Code :**  
LBNY-00017-11000280



**Description :**

Constructed of heavyweight spunbonded polypropylene  
Designed to provide strength and fluid-resistance while offering breathability and comfort.  
Synthetic elastic conforms to head for a personalized fit  
24 in. diameter  
100% latex-free.

**Technical Specification :**

Blue Premium Bouffant Cap



**Laboratory Instrument India**