



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Black Series Heavy Duty Vest

Product Code :
LBNY-00017-11000281



Description :

Features zipper front closure reinforced with durable webbing. Premium Black Series Heavy Duty Vest is constructed of 100% polyester mesh.
Reflective material (2 in. wide) with 3 in. breathable contrasting color
Left and right chest mic tab
Ultra-Cool 100% polyester mesh construction
Zipper front closure
Reflective trim around arms and piping down sides for added low light visibility (Except for size: Medium)
Pockets:
Left chest heavy duty gusseted radio/utility pocket with reflective piping
Two inside patch pockets
Two lower outside cargo pockets with adjustable flaps and grommets
Right chest 2-tier pencil pocket.

Technical Specification :

Premium Black Series Heavy Duty Vest



Laboratory Instrument India

