



Email : sales@laboratoryinstrumentindia.com

Product Name :
Unisex Lab Coats with Knit Cuffs

Product Code :
LBNY-00017-11000282



Description :

Made of 65/35 polycotton
One outside left breast pocket and two outside lower pockets
Standard pointed collar
The knit wrist cuffs of these Unisex Lab Coats provide a comfortable and secure fit.
Knit cuffs
Front snap closure
White.

Technical Specification :

Unisex Lab Coats with Knit Cuffs



Laboratory Instrument India