



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Unisex Lab Coats with Knit Cuffs

**Product Code :**  
LBNY-00017-11000282



**Description :**

Made of 65/35 polycotton  
One outside left breast pocket and two outside lower pockets  
Standard pointed collar  
The knit wrist cuffs of these Unisex Lab Coats provide a comfortable and secure fit.  
Knit cuffs  
Front snap closure  
White.

**Technical Specification :**

Unisex Lab Coats with Knit Cuffs



**Laboratory Instrument India**