



Email : sales@laboratoryinstrumentindia.com

Product Name :
Unisex Specialized Cuffed Lab Coat

Product Code :
LBNY-00017-11000283



Description :

This lab coat is off the charts with all our best material production features. Featuring a 6-gripper closure, neck gripper, and a chest pocket, this Unisex Specialized Cuffed Lab Coat looks as professional as it is durable. Lab, medical facility, and clinic workers will appreciate the knit cuffs that keep sleeves from interfering with work and extended back coverage to protect from messy jobs. Plus, who doesn't need the storage offered by a chest pocket and two lower pockets? Put this lab coat to the test and see why it's a customer favorite.

Technical Specification :

Unisex Specialized Cuffed Lab Coat



Laboratory Instrument India