



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cleanroom Badge Clips

**Product Code :**  
LBNY-00017-11000285



**Description :**

Ideal for visually identifying different divisions, departments, titles, duties, or access categories.  
Double-bagged to "Class 10" specifications, for use in critical environments  
Secure badges or ID cards to garments. Cleanroom Badge Clips are made of three inch poly for badges or garments with no pockets. These swivel clips are metal-free.  
Non-metal clips perfect for metal-free environments.  
Contain no metal  
50/pack.

**Technical Specification :**

Cleanroom Badge Clips



**Laboratory Instrument India**