



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Laboratory Scale Spray Dryer

**Product Code :**  
LBNY-0001-990010



**Description :**

Laboratory Scale Spray Dryer

**Technical Specification :**

Individual controller for temperature, vacuum and feed rate using control knobs.

Spray dry with aqueous as well as solvent based solutions.

Compact table top model; enables spray drying trial with even 30 ml samples.

Designed with a SS material & reinforced with flameproof motors, to withstand the effect of acid & solvent spray drying.

Simple to operate and clean, economical and gives quick results.



**Laboratory Instrument India**