



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bouffant Cap and Beard Cover Hood

Product Code :
LBNY-00017-11000286



Description :

Economical, breathable spun bound polypropylene fabric
Elastic around the facial opening and neck ensures a secure fit
This is manufactured from a specially formulated spunbond polypropylene (SBPP) fabric, and can be used in cleanroom environments.
One size fits most.

Technical Specification :

Bouffant Cap and Beard Cover Hood



Laboratory Instrument India