



Email : sales@laboratoryinstrumentindia.com

Product Name :
QC Sleeved Apron

Product Code :
LBNY-00017-11000290



Description :

Knee-length apron/gown design with attached sleeves
Lightweight and durable fabric
Provides enhanced frontal protection
Tightly sewn seam is covered with garment fabric to reinforce seam and reduce potential for particle penetration
Open back for wearer comfort
Bound seams around neck opening with snap closure
Elastic opening for tighter fit at wrist
QC fabric provides at least 30 minutes of protection against 42 chemical challenges. It is used for light splash protection in a variety of industrial environments, including petroleum refining, pulp and paper manufacturing, food processing, chemical processing, and pharmaceutical manufacturing.
Ties at waist and elastic wrists
Yellow colored apron is 44 in. long for neck-to-knee coverage.

Technical Specification :

QC Sleeved Apron



Laboratory Instrument India

