



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Gyratory Shaker

**Product Code :**  
LBNY-0001-1000001



#### Description :

Gyratory Shaker, Gyratory Shaker Equipments, Gyratory Shaker Tools, Gyratory Shaker Tool Kits, Gyratory Shaker Manufacturers, Gyratory Shaker Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara

#### Technical Specification :

Gyratory Shaker 5-70RPM 3-12° THROW 230V 50/60HZ.

Sieve Shaker Gyratory is designed to carry up to 7 sieves of 6 inch or 8 inch diameter. The 12 inch diameter sieves can be used with the help of adaptor. A 1/4 HP. motor through a reduction gear immersed in oil drives the shaker.

The sieve table does not rotate but is inclined from the vertical axis and the direction of inclination changes progressively in clockwise direction. If the stop-pin below the table is removed, the shaker can have a rotary motion. In addition to this gyratory motion of the table, there is an upward and downward movement ensuring that each square cm of the sieve is utilized.

A pair of rods and a holder is supplied. The holder can be fixed on the top of the upper most sieves, and thus the sieve set is firmly held. Suitable for operation on 230 volts, 50 cycles single phase A.C Supply. Adaptor for 30cm dia. test sieves can be fitted in above.



**Laboratory Instrument India**