



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Polynit Heat Seal Wipes

**Product Code :**  
LBNY-00017-11000296



**Description :**

Polynit Heat seal Wipes are made of 100% knitted textured polyester for very low levels of particles and fibers. Polyester wipes are chemically resistant and exceptionally low in particles and fibers. Interlock knit creates a durable fabric, the addition of a periodic additional stitch of a no-run interlock knit prevents the fabric unraveling creating an even more durable fabric with strong edges. Good sorbency with solvents and are abrasion and chemical resistant. 100% laundered knitted polyester creates very low levels of particle and fibers. Chemically resistant and exceptionally low in particles and fibers.

**Technical Specification :**

Polynit Heat Seal Wipes



**Laboratory Instrument India**