



Email : sales@laboratoryinstrumentindia.com

Product Name :
Non-Sterile Integrated Hood and Mask

Product Code :
LBNY-00017-11000298



Description :

Hood and mask combination guards against the risk of contamination in cleanroom environments.
Stretch fit elastic headband for more secure fit
Pull-down ties in back to assist in the gowning process
Integrated hood and mask sonically welded together combines two donning steps into one
Universal size
White
Latex- and rubber-free
Bulk packed.

Technical Specification :

Non-Sterile Integrated Hood and Mask



Laboratory Instrument India