



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Laboratory Constant-Temperature Drying Oven

**Product Code :**  
LBNY-0001-1000002



**Description :**

Laboratory Constant-Temperature Drying Oven

**Technical Specification :**

Laboratory Constant-Temperature Drying Oven 70L capacity  
External material: cold-rolled steel with anti-bacteria powder coating.  
Stainless steel inner chamber, round angle structure, adjustable shelves.  
Double-layer glass observation window, adjustable airtight buckle lock.  
PID control with LED display  
Over-temperature protection.



**Laboratory Instrument India**