



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Constant-Temperature Drying Oven

Product Code :
LBNY-0001-1000002



Description :

Laboratory Constant-Temperature Drying Oven

Technical Specification :

Laboratory Constant-Temperature Drying Oven 70L capacity
External material: cold-rolled steel with anti-bacteria powder coating.
Stainless steel inner chamber, round angle structure, adjustable shelves.
Double-layer glass observation window, adjustable airtight buckle lock.
PID control with LED display
Over-temperature protection.



Laboratory Instrument India