



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Clean Full Face Opening Hoods

**Product Code :**  
LBNY-00017-11000313



**Description :**

Clean open face hoods provide an inherent barrier to particles, microorganisms, and non-hazardous light liquid splash. Available in non-sterile, clean-processed, and sterile versions.

Full face opening

Ties with loops for adjustable fit

Bound seams are covered with garment fabric for reinforcement

Bound hood opening for lower particle shedding

One size fits most

Clean-processed garments Cleanroom 4 standards

Available in bulk packaging with double transparent poly liners

Available in non-sterile, clean-processed, and sterile versions.

**Technical Specification :**

Clean Full Face Opening Hoods



**Laboratory Instrument India**

