



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vinyl Closed-Back Aprons, Plain edge

Product Code :
LBNY-00017-11000314



Description :

With heat-sealed elastic cuffs
6mil thick aprons have electronically welded binding strip around neck edge only; other edges are unhemmed
8mil thick aprons have electronically welded vinyl binding strip that covers all edges for smoother finish
Available with snap -on nylon ties
Neck ties, 15 in. long (38cm); waist ties, 29 in. long (74cm)
Blue.

Technical Specification :

Vinyl Closed-Back Aprons, Plain edge



Laboratory Instrument India