



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Tube Knit Sleeves with Aramid Fiber

**Product Code :**  
LBNY-00017-11000315



**Description :**

Made from 100% Aramid fibers  
Feature thumb slot to hold sleeve down for a more comfortable fit  
Provide workers' palms and arms with the same cut and heat protection as aramid fiber string knit gloves  
Launderable.

**Technical Specification :**

Tube Knit Sleeves with Aramid Fiber



**Laboratory Instrument India**