




Email : sales@laboratoryinstrumentindia.com

Product Name :
Responder CSM Suits

Product Code :
LBNY-00017-11000317



Description :

Hood features a large 40mil PVC face shield
Expanded back allows room for SCBA
Two exhaust valves reduce internal pressure buildup
Replaceable butyl gloves attached with ring-and-clamp system
Triple overlapping storm flap
Level A Fully Encapsulating (Gastight) Suit
Front-entry suit with gastight PVC zipper provides maximum protection for skin and eyes
Attached sock boots with boot flap
Double-sealed seams and 5mil PTFE face shield overlay
Limited-use suits are designed specifically for dealing with deadly chemical-warfare agents, or Chemical Surety Materials (CSM).
Responder CSM fabric is made to military specifications, with 100% inflation testing plus lot testing of suit construction and fabric.
Butyl gloves attached with jam-fit glove system
Attached sock boots with splash guards
Double-sealed seams
Extended double-storm flap works with hood to accommodate respirator mask
Weight: 9 lb. (4.1kg)
Level B Fully Encapsulating (Non-Gastight) Coverall
Zipper front with attached respirator hood
Can be used with optional over hood to protect SCBA
Weight: 5 lb. (2.3kg).

Technical Specification :

Responder CSM Suits



Laboratory Instrument India