




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Responder CSM Suits

**Product Code :**  
LBNY-00017-11000317



**Description :**

Hood features a large 40mil PVC face shield  
Expanded back allows room for SCBA  
Two exhaust valves reduce internal pressure buildup  
Replaceable butyl gloves attached with ring-and-clamp system  
Triple overlapping storm flap  
Level A Fully Encapsulating (Gastight) Suit  
Front-entry suit with gastight PVC zipper provides maximum protection for skin and eyes  
Attached sock boots with boot flap  
Double-sealed seams and 5mil PTFE face shield overlay  
Limited-use suits are designed specifically for dealing with deadly chemical-warfare agents, or Chemical Surety Materials (CSM).  
Responder CSM fabric is made to military specifications, with 100% inflation testing plus lot testing of suit construction and fabric.  
Butyl gloves attached with jam-fit glove system  
Attached sock boots with splash guards  
Double-sealed seams  
Extended double-storm flap works with hood to accommodate respirator mask  
Weight: 9 lb. (4.1kg)  
Level B Fully Encapsulating (Non-Gastight) Coverall  
Zipper front with attached respirator hood  
Can be used with optional over hood to protect SCBA  
Weight: 5 lb. (2.3kg).

---

**Technical Specification :**

Responder CSM Suits



**Laboratory Instrument India**