



Email : sales@laboratoryinstrumentindia.com

Product Name :
Clean and Sterile Bouffant

Product Code :
LBNY-00017-11000319



Description :

Bouffant creates a barrier to dry particles, microorganisms, and nonhazardous liquids. Suitable for ISO 5 (class 100) controlled environment.

Elastic headband

Comfortable, lightweight and durable

Garments available gamma sterilized to an SAL of 10⁻⁶

Serged seams

Color: White.

Technical Specification :

Clean and Sterile Bouffant



Laboratory Instrument India