



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Protective Coveralls: With Collar

**Product Code :**  
LBNY-00017-11000321



**Description :**

Breathable, cloth-like material treated to provide antistatic properties for use in controlled environments  
Proprietary bilaminate composite fabric provides superior liquid penetration  
Compared favorably to leading fabric in physical property tests covering resistance to ripping and tearing.  
Affordable high-barrier protection  
Lightweight coveralls (1.8 oz. per sq. yd.) provide excellent barrier protection against nontoxic particle, liquid and spray hazards. True aerosol- and liquid-proof microporous barrier demonstrates >99% holdout of particles down to 0.35 $\mu$ m.  
Long torso with polyester front zipper measures 26 in. long (66cm)  
White.  
Low linting for use in critical applications.

**Technical Specification :**

Protective Coveralls: With Collar



**Laboratory Instrument India**

