



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pro Clean Boot Covers

Product Code :
LBNY-00017-11000323



Description :

Serged seams have multiple interlocking threads that are sewn around the raw edges of garment material to create a strong, stress-resistant seam
Boot covers are 18 in. with covered elastic top and embossed PVC soles.
High-top boot covers with elastic tops, and ankles
Height: 18 in. (45.7 cm)
Color: white.

Technical Specification :

Pro Clean Boot Covers



Laboratory Instrument India