



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chemical Protective Apron with Sleeves

Product Code :
LBNY-00017-11000325



Description :

Chemical Protective Apron with Sleeves for exceptional chemical protection
Double cuff with knitted inner cuff
Welded and taped seams
Velcro fastening to neck
Tie fastening at the waist
Applications:
Contaminated environments
Tranferring liquids and solids
Opening and draining pumps valves or lines
Loading and unloading of process equipment
Chemical handling
Filling, blending and charging of raw materials
Hazmat response
Assembly and disassembly of large pipes.
Raw material sample collection
Concentration mixing and palletisation.

Technical Specification :

Chemical Protective Apron with Sleeves



Laboratory Instrument India