



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Chemical Protective Apron with Sleeves

**Product Code :**  
LBNY-00017-11000325



**Description :**

Chemical Protective Apron with Sleeves for exceptional chemical protection  
Double cuff with knitted inner cuff  
Welded and taped seams  
Velcro fastening to neck  
Tie fastening at the waist  
Applications:  
Contaminated environments  
Tranferring liquids and solids  
Opening and draining pumps valves or lines  
Loading and unloading of process equipment  
Chemical handling  
Filling, blending and charging of raw materials  
Hazmat response  
Assembly and disassembly of large pipes.  
Raw material sample collection  
Concentration mixing and palletisation.

**Technical Specification :**

Chemical Protective Apron with Sleeves



**Laboratory Instrument India**