



Email : sales@laboratoryinstrumentindia.com

Product Name :
SL Apron

Product Code :
LBNY-00017-11000328



Description :

Made of spunbonded olefin laminated with provide chemical barrier
Provides enhanced frontal protection
Knee-length apron/gown design with attached sleeves
Open back for wearer comfort
Bound seams around neck opening with snap closure
Ties at waist
Taped seams provide strong chemical resistance against heavy liquid splashes
Sewn seam is covered with strip of compatible chemical-resistant material through heat sealing
Elastic opening for tighter fit at wrist
Available in multiple sizes
44 in. long for neck-to-knee coverage
Color: White.

Technical Specification :

SL Apron



Laboratory Instrument India

