



Email : sales@laboratoryinstrumentindia.com

Product Name :
Endurosaf Aprons

Product Code :
LBNY-00017-11000332



Description :

Strong, tough, 8 to 9 mil die cut polyurethane material offers a unique combination of characteristics to outperform neoprene, nitrile, vinyl and other films
High performance, adaptable protection for a wide variety of jobs and are resistant to abrasion, chemicals, oils and fats
Soft, comfortable, easy to clean and remain flexible in cold temperatures.
Sealed ties and adjustable neck straps provide a durable, consistent fit for a wide variety of users.

Technical Specification :

Endurosaf Aprons



Laboratory Instrument India