



Email : sales@laboratoryinstrumentindia.com

Product Name :
A20 Lab Coats

Product Code :
LBNY-00017-11000335



Description :

Snap front closure with snap at collar
Long length
Elastic wrists
Professional A20 Lab Coats feature breathable, patented MICROFORCE barrier SMS fabric
No pockets
Typical uses: Laboratories; Decontamination modesty garment; Remediation; Pesticides; Heavy dust environments.
Passes NFPA 99 criteria for antistatic materials.

Technical Specification :

A20 Lab Coats



Laboratory Instrument India