



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Autoclave 18L capacity

Product Code :
LBNY-0001-1000005



Description :

Laboratory Autoclave 18L capacity

Technical Specification :

Easy operation, safe and reliable.
The sealing ring is made of silicon rubber.
Two safety valves keep double safety (for timing control type).
It has self-adjustment function (for timing control type)
It is well sealed and effectively prevents steam leakage.
Over temperature & over pressure auto-protection.
Electric heating or gas heating can be used (for Dual fuel Type).
The sterilizer is provided with a power switch, user can control the operation easily.
Immersion type electric heating pipe, thermal efficiency is high, stainless steel material tube is durable to use.
Steam-water inner circulation system: No steam discharge and the environment for sterilizing will be clean and dry.
Double groove pressure gauge, could clearly indicate, convenient observation, reliable performance.



Laboratory Instrument India

