



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Laboratory Autoclave 18L capacity

**Product Code :**  
LBNY-0001-1000005



**Description :**

Laboratory Autoclave 18L capacity

**Technical Specification :**

Easy operation, safe and reliable.  
The sealing ring is made of silicon rubber.  
Two safety valves keep double safety (for timing control type).  
It has self-adjustment function (for timing control type)  
It is well sealed and effectively prevents steam leakage.  
Over temperature & over pressure auto-protection.  
Electric heating or gas heating can be used (for Dual fuel Type).  
The sterilizer is provided with a power switch, user can control the operation easily.  
Immersion type electric heating pipe, thermal efficiency is high, stainless steel material tube is durable to use.  
Steam-water inner circulation system: No steam discharge and the environment for sterilizing will be clean and dry.  
Double groove pressure gauge, could clearly indicate, convenient observation, reliable performance.



**Laboratory Instrument India**

