



Email : sales@laboratoryinstrumentindia.com

Product Name :
10 Basic Disposable Lab Coats

Product Code :
LBNY-00017-11000336



Description :

Designed for workers in a range of industries for nonhazardous dry particle applications.
Spunbond meltblown spunbond (SMS) material
Three pockets (one left chest pencil, two lower front)
Disposable lab coats that are optimized for comfort, softness and breathability
Set sleeves
Snap closure (six plus one adjustable)
Serged seams
Knit collar
Knit cuffs.

Technical Specification :

10 Basic Disposable Lab Coats



Laboratory Instrument India