



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
F Taped Seam Coveralls

**Product Code :**  
LBNY-00017-11000338



**Description :**

F coverall is made of barrier film laminated to substrate  
For use when potential exposure to industrial chemicals and chemical warfare agents exists  
Widely used by military personnel and first responders for chemical warfare situations  
Provides protection against over 180 chemical challenges for at least 30 minutes  
Strong and durable coverall with broad chemical barrier  
Front zipper closure and storm flap with tape closure  
Elastic wrists and elastic ankles.  
Available in low-visibility gray fabric for discretionary purposes  
With taped seams and attached hood (respirator fit).

**Technical Specification :**

F Taped Seam Coveralls



**Laboratory Instrument India**