



Email : sales@laboratoryinstrumentindia.com

Product Name :
Aramid Fiber Assembler 100 Sleeve

Product Code :
LBNY-00017-11000339



Description :

Features knit wrist thumb slot to protect palm and wrist
With thumb slot for extra protection at palm and wrist
Flame-retardant for use in areas with moderate heat
100% Aramid fibers interlock knitted fabric for superior protection
Offers moderate cut and abrasion protection for light-duty applications
Sleeves provide cut and abrasion resistance suitable for most general assembly applications.
Length: 22 in. (56cm).

Technical Specification :

Aramid Fiber Assembler 100 Sleeve



Laboratory Instrument India