



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Centrifuge 15ml x 8 Bucket

Product Code :
LBNY-0001-1000006



Description :

Laboratory Centrifuge 15ml x 8 Bucket

Technical Specification :

Economical unit for small capacity centrifuging requirements.
The standard unit is fitted with resiliently mounted motor with 3500 RPM for vibration free performance, a built in 5 step speed regulator and one non removable angle motor.
Complete with tube carriers.



Laboratory Instrument India