



Email : sales@laboratoryinstrumentindia.com

Product Name :
ESD Cleanroom Overshoes

Product Code :
LBNY-00017-11000347



Description :

Electrostatic Discharge (ESD) cleanroom overshoes
Constructed of non-woven spunbonded polypropylene
ESD cleanroom overshoes have non-marking conductive tape, providing dissipative properties when worn as instructed.
Thickness 0.08 mm
Choose white 16 in. or 18 in. or blue 18 in.

Technical Specification :

ESD Cleanroom Overshoes



Laboratory Instrument India