



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
224 Series Lab Coats

**Product Code :**  
LBNY-00017-11000353



**Description :**

Light and durable lab coats provide upper body coverage, while staying comfortable.

Protect against particulates and other contaminants in controlled environments.

Lowest level particle shedding

Serged seams

Five-snap closure

Spinning process creates barrier to particles, micro-organisms and nonhazardous liquids

Laydown collar

Raglan sleeves

Two lower front pockets and one chest pocket.

Applications:

Use in controlled environments with high standards for particle and microbiological contamination control;

Pharmaceutical, biotech, medical device, electronic settings.

**Technical Specification :**

224 Series Lab Coats



**Laboratory Instrument India**

